

Steps for success!

- Alcohol-based products may interfere with testing. These products include but are not limited to the following: alcohol-based mouthwash, cough syrup, hand sanitizer, colognes, perfumes, hair spray, or similar products. This is the most common cause of contaminated results. Always wait at least 5 minutes after drinking, eating, or smoking to ensure the most accurate results.
- Do not blow smoke, saliva, or other contaminants into the mouthpiece, as the sensor may become damaged, requiring a replacement at your expense.
- Although the i70 pro breathalyzer has been rigorously tested for evidential use, it is essential to understand that it is only a screening device and should be treated as such in this application. Exercise caution, and if you have consumed any amount of alcohol, never drink and drive.
- Keep the i70 pro device out of extreme temperatures and out of the reach of children.
- Please use alkaline batteries only.
- To maintain the accuracy of the i70 pro, re-calibrate your device every 1000 tests or at least every 12 months.. To recalibrate, you must contact the organization you purchased your device from for details and related fees.

Troubleshooting Tips

Connectivity Issue:

Hand Held Breathalyzer won't connect to the app

When the app and the device are properly connected via Bluetooth, the i70 pro will display "APP" on the breathalyzer, and the app will display "Start Test" in the status window.

If this doesn't occur, complete the following steps:

- Power off the breathalyzer and close the app. Then reopen the app, and power on the breathalyzer
- If the problem persists, with the App open, toggle the phone's Bluetooth off, pause for about 15 seconds, power on the i70 Pro, then turn Bluetooth back on
- If the problem persists, make sure the phone is not connected to another Bluetooth device, which could prevent connection
- If the problem persists, ensure all app permissions are "allowed" on the smartphone
- If the problem persists, ensure Location settings are set to "always allow" or "allow while using the app"
- If the problem persists, ensure there are no other phones with applications nearby (if so, temporarily turn off their Bluetooth)

If any of the above are true, correct as instructed, force the app closed, reopen the app, and power up the handheld breathalyzer.

- If the problem persists, confirm that the battery in the Handheld Breathalyzer is sufficient, and replace the battery if necessary
- If the problem persists, Toggle from Wi-Fi to Cellular or Vice Versa (depending on the current form of data transmission)
- If the problem persists, log out and then log back into the app
- If the problem persists, uninstall the app, then reinstall, and log back in
- If the problem persists, schedule breathalyzer for Calibration (inspection)

How to replace batteries:



1. Use a penny or flathead screwdriver at the device's bottom notch
2. Open the battery cover slightly from the bottom to the upper side by pulling the battery cover from the unit
3. Remove the battery cover completely, then replace the batteries. Re-assemble the battery cover in reverse order of disassembly.

Error Messages

Message	Instruction
Power off i70 pro	Hold the power button down for approximately 5 seconds, the unit will display "off" with a countdown 3, 2, 1, 0, then power off.
(BAT): Low Battery Indicator	This can occur in the in app, or the Hand Held Breathalyzer. Confirm the battery in the handheld breathalyzer is fresh, if necessary replace the battery
(FLOW): Blow Error	The device detected an issue while you were presenting your sample. This could be from blowing too softly, not continuously, not long enough, or incorrectly. Power off the handheld breathalyzer, wait a few seconds, power back on, then try another sample
(ERR1) : Pressure sensor error –	Ensure you are not providing sample too early, and you are blowing firmly, steady and consistently
(ERR2) : Sensor error –	Contact your Service Provider
(ERR3) : Temperature Error (Out of Operational temperature range)	This indicates that you are using your breathalyzer in extreme temperature conditions. Take your breathalyzer to a climate controlled area and allow it to warm up or cool down. Once the unit temperature is between 23° F and 120° F return to the vehicle and attempt to test. To avoid this error always take your breathalyzer with you when not in your vehicle.
(CAL) : Calibration Error	The device is required to be calibrated on regular intervals (the shorter between 12 months, or 1000 tests) if you receive this error, you have exceeded one of these thresholds, and will need to have your device calibrated iSOBER 70 will display "CAL" sign for 3 seconds then warming up for test.

If any problem persists, the device must be returned for calibration/inspection. Contact service provider to make those arrangements. Replacement device arrangements should be made based on warranty status.